



3-QUINUCLIDINYL BENZILATE (BZ)

Agent Overview:	3-Quinuclidinyl Benzilate (BZ) is considered an incapacitating and hallucinogenic agent. The onset of incapacitation is dose-dependent; it might occur as early as 1 hour after exposure and continue up to 48 hours. BZ poisoning can occur by inhalation, ingestion or skin absorption. BZ poisoning results in an anticholinergic syndrome consisting of a combination of signs and symptoms that might include hallucinations, agitation, dilated pupils (mydriasis), blurred vision, dry, flushed skin, urinary retention, ileus, tachycardia, hypertension, elevated temperature ($>101^{\circ}\text{F}$), restlessness, dizziness or giddiness, failure to obey orders, erratic behavior, stumbling or staggering, vomiting, disrobing, picking behavior, mumbling, slurred or non-sensical speech, stupor, or coma.
Lab Specimen Submittal:	Authorization by Dr. Jane Getchell, Director DPHL, Christina Pleasanton, Deputy Director DPHL, or, Tara Lydick, Chemical Terrorism Coordinator
Type of Specimens Necessary for Testing:	<u>Urine:</u> 25mL-50mL screw-cap container <u>Whole blood:</u> <ul style="list-style-type: none">- Three, 5 or 7 mL purple-top (EDTA) tubes (vacuum fill only) AND- One, 5 or 7 mL gray-top or green-top (vacuum fill only) AND- Two empty urine cups and two empty blood tubes for each lot # submitted.
Packaging Instructions:	<u>Urine:</u> freeze specimen upon collection and ship on dry ice; label specimens and package individually; wrap with absorbent material. <u>Whole blood:</u> place filled tubes in a shipping rack for stability, wrap the tubes in absorbent material, and secure with tape. Place the packaged tubes in a secondary unbreakable container. Specimens may be shipped cold but <i>must not</i> be frozen. <u>All:</u> wrap sample containers and layers with tamper-proof, waterproof forensic evidence tape initialed by collector, being careful not to cover the sample labels. Ship to DPHL as Biological Substance, Category B (Packing Instruction 650).
Receiving Specimens at DPHL:	Authorized specimens are received at the lab's rear (loading dock) entrance. Submitters must present ID and complete proper forms.
Handling of Specimens:	Chain of custody procedures must be followed. Urine specimens should be frozen upon collection.
Reporting Results:	Reference laboratory will notify DPHL as soon as possible. The Potential Exposure Call Log and Information form will provide information for notification of testing results outside of the DPHL.
Contact Information:	Jane Getchell, DrPH, Director DPHL, 302-223-1520 (24/7 coverage) 24/7 Emergency Contact Number: 1-888-295-5156 Revised: 04/08