



Division of Public Health




Delaware's Public Health Laboratory and Sentinel Laboratories' Partnership

**Marion Fowler, MT(ASCP)
Microbiologist II**

Lab Preparedness



- CDC/ LRN Preparedness Grant 2001
- Public Health Labs – Increase capacity to confirm BT agents/molecular tests
- Sentinel Lab – need Rule/Out training
- Packaging & Shipping infectious agents
- Laboratory Preparedness Advisory Committee (LPAC)



Proficiency testing


- A means to test sentinel labs on their ability to rule/out BT agents
- Most PHL's are not in PT business
- CAP (College of American Pathologists) already sends PTs to most labs
- LPS survey was not testing appropriate organisms or stains. Answer choices not applicable
- APHL interceded to assist in getting CAP to work with PHL staff on revising format.

LPX CAP Survey



- CAP & PHL workgroup revised survey
- List of surrogate organisms
- List of appropriate responses
- Contact LRN reference lab for r/o testing
- Packaging & shipping exercise

DE CAP LPX-A April 2011



- LPX-01 (*non-BT culture*)
- 10 reporting labs- 9 sentinel, 1 reference
 - 10 labs: non-BT culture
 - Gram negative rod, normal growth at 24 hrs for a gnr, growth on SBA, CHOC and MAC, oxidase pos, indole pos, urea neg, growth at 42°, arginine neg

—


Delaware Results

- LPX-02 (*Bacillus anthracis*)
 - 1 lab DPHL: confirmed *Bacillus anthracis*
 - 9 labs: unable to rule out *Bacillus anthracis* using ASM Sentinel Laboratory Guidelines:
 - non-motile (motility/TTC)
 - non-hemolytic
 - catalase positive
 - large gram positive rods



Delaware Results

- LPX-03 (*Francisella tularensis*)
 - 1 lab DPHL: confirmed *Francisella tularensis*
 - 9 labs: unable to rule out *Francisella tularensis* using ASM Sentinel Laboratory Guidelines:
 - oxidase negative
 - urea negative
 - catalase- slightly positive
 - beta-lactamase positive
 - slow-growing, extremely tiny faintly staining gncb



Acceptable Responses and Contact

- Reported for LPX-A 2011
 - LPX-01 DE: 100% reported non-BT culture
 - LPX-02 DE: 100% reported unable to rule out *Bacillus anthracis*
 - LPX-03 DE: 100% reported unable to rule out *Francisella tularensis*
- Contact with DPHL
 - 2011 LPX-A: 100% DE sentinel labs
 - 2010: LPX-A and LPX-B: 100% DE sentinel labs
 - 2009: LPX-B: 78% of DE sentinel labs

Packaging & Shipping

- All labs contacted DPHL by the due date of the survey.
- All labs received on-site training in the fall of 2010; therefore, no labs were required to package and ship at this time.
- All labs will be required to package and ship for the LPX-B survey in the fall of 2011.



Agents of Bioterrorism: Annual Sentinel Lab Update



- 2002: First Agents of Bioterrorism Workshop was held in 2002
- Updates continued yearly until!!!!
- 2006: Train-the-trainer NLTN course available which included hands-on training
- 2007,2008 and 2010: training available which included the popular “wet workshop”

Full Day Training includes:



- Review of Sentinel Laboratory Guidelines for Suspected Agents of Bioterrorism (ASM, CDC and APHL)
- DPHL protocols and forms
- Safety training
- “Wet Workshop” in afternoon using vaccine strain or attenuated bioterrorism organisms. Each BT organism, along with two other organisms that closely mimic a BT organism are plated on a variety of agars and plates are available to read at 24, 48 and 72 hours of incubation.

Conclusions

- BT workshop training is necessary at least every other year due to high turnover in sentinel laboratories
- Best way to learn is through a combination of wet workshops, review of Sentinel Laboratory Guidelines for Suspected Agents of Bioterrorism and participation in the CAP LPX surveys

Conclusions

- BSL2/3 training given during the workshop can reduce the risk of an acquired laboratory infection when working with organisms transmitted by aerosols
- Packaging & shipping training is necessary for several reasons:
 1. Sentinel laboratories must be able to properly package an isolate for rule out of a BT agent. This is for the safety of the courier and the laboratory personnel at both laboratories.
 2. In case of a public health emergency, the sentinel laboratory may be requested by CDC or DPHL to send cultures/specimens, etc. directly to CDC.

Questions?



marion.fowler@state.de.us or debra.rutledge@state.de.us

**Delaware Public Health Laboratory
30 Sunnyside Road
Smyrna, Delaware 19977
302- 223-1520**